



Glass Pipette Preparation

Ethan Bai




Introduction

Many methods of injecting therapeutics to the brain, such as

- Intracerebral injections

Has its upsides, such as its precise targeting, however one major negative is the its invasive nature



Motivations

Thin tipped tools



- Thin-tip glass pipette

Problems

- Affordability
 - Price of obtaining glass pipettes
 - Requires conversion of professional programmable pipette pullers (thousands of dollars)
- Accessibility
 - Non electrophysiological labs lack access to pipette pullers

3D Printing can drastically reduce the price

Pulling Device




My Project: Processing and Calibrating Glass Pipettes

- Calibrating
- Trimming
- Beveling

Calibrating

- Different therapeutics to be injected are different sizes, so we must cater our pipettes to these needs
 - Liquid delivery - 20 μm
 - Cell delivery - 40-50 μm
- Determining where to cut pipette
 - Visible measuring scale in the microscope lens
 - Reticle
 - Allows measurement of items visible in the microscope



Reticle



Reticle

Microsoft Word used to create the reticle

1. Measure the diameter of the microscope lens
2. Create a circle of that diameter in word
 - Set height and width equal to the diameter
3. Convert desired length in microscope to actual length on the reticle
 - For 100x magnification, $20\text{ }\mu\text{m} * 100 = 2\text{ mm}$
 - $40\text{ }\mu\text{m} * 100 = 4\text{ mm}$
4. Create line of that length and center it in the circle with align tool
 - Multiple lines -> group first
5. Print out reticle on a clear film
6. Cut out the reticle
 - Cut out "ears" on the reticle so it's easy to take out of the microscope



Problem with Old Reticle

Problem:

- Intended lengths obtained from calculations were off
 - Ex: The line for 20 μm was actually 150 μm through the microscope
 - Hidden lens within the microscope

Solution:



1. Calculate magnification factor
 - Ex: $150/20 = 7.5$
2. Divide length of the line on the reticle by magnification
 - $2 \text{ mm}/7.5 = 0.2667 \text{ mm}$

Trimming

Shining a bright light through one of the lenses projects a bright dot

Trimming Procedure:

1. Match the line on the reticle with the pipette to measure desired diameter
2. Turn off the microscope light and turn on the flashlight
3. Make a cut on the pipette where the bright dot is



Beveling

We can bevel the pipettes with a spinning hard drive disk

Put the tip on the surface of the disk very carefully at a small angle



- Too much force could break the tip
- Too large of an angle could result in a dull tip



Software

We can confirm if our measurements were correct using an app called imageJ.

1. Take screenshot of the grid on the slide
 - Each space is $10\text{ }\mu\text{m}$
2. Draw a line from origin to any line on the scale
3. Calculate that distance in μm
4. Analyze -> Set Scale
 - Change known distance to calculated length
 - Change units to μm
 - Check global



Thank You
Questions?